



Glass ionomer cements & lining and base materials

Jayden Yoon, DDS

GIC Intro


- Consists of glass particles that integrate into a network of acidic polymers through acid-base reaction
- Used in most disciplines in dentistry
- Benefits: fluoride releasing, natural adhesion (chelation) to enamel and dentin




GIC components

- Glass
 - Alumina (Al₂O₃) and silica (SiO₂) with fluoride salts
 - Impacts ion release and esthetics
 - Powder form incorporates larger particles for filling GIC vs. smaller for lining GIC
- Polymeric acid (polyacid)
 - Influence setting strength: higher MW --> stronger set
- Water
- Reaction modifiers
 - Tartaric acid most common


Acrylic acid




Itaconic acid




Maleic acid



Tartaric acid



Phosphonic acid



GIC forms

- Hand-mixed powder-liquid: strict ratios, mix timing sensitive
- Encapsulated: activation by breaking a membrane separating powder/liquid, then triturating





TABLE 2.3.1

Comparison of Properties of Various Formulations and Versions of Glass Ionomers

Property (Unit)	Hand-Mixed GIC	Encapsulated GIC	Hand-Mixed RMGIC	Encapsulated RMGIC
Microhardness (Knoop hardness number) at 1 day	40	58	33	49
Microhardness (Knoop hardness number) at 30 days	47	63	28	37
Biaxial flexural strength (MPa) at 1 day	42	56	123	136
Biaxial flexural strength (MPa) at 30 days	48	52	92	174
Compressive strength (MPa) at 1 day	159	205	165	170
Compressive strength (MPa) at 30 days	156	193	171	182


GIC setting reaction

- Acid-base reaction between polymeric acid and the glass
- Calcium and aluminum ions are released from glass to form a salt matrix
- Sodium fluoride is released and flows throughout restoration
- 3 phases
 - Dissolution
 - Gelation
 - Forming




Setting reaction: dissolution

- Acid reacts with outer layer of glass
 - Removes metal ions (aluminum, silica, sodium, fluorine, calcium)
 - Only silica gel remains
 - Hydrogen ions from carboxyl groups diffuse to glass




Setting reaction: gelation

- Calcium ions chelate polymeric chains rapidly
 - Aluminum chelates at a slower pace due to trivalency
- Critical notes
 - Needs water to react (do not overly dry restoration)
 - Too much water weakens material and dulls esthetics (reduce water contamination)
 - Aluminum may diffuse out in liquid (cannot form chains)




Setting reaction: hardening

- Aluminum leads the hardening reaction with its trivalency (more cross-linking), but is slower than Ca
 - ~24 hrs till completion: hardness, strength, translucency, and esthetics improve
 - Reaction believed to continue for weeks, some dentists add on a layer of varnish to supplement




Tartaric acid and water effects

- Tartaric acid
 - Extends working time by slowing down release of Ca ions
 - Then speeds up reaction by increasing rate of reaction of Al cross-linking
- Water
 - Allows ions to migrate, complete reaction, and fluoride movement



Adhesive properties

- Adhesion to tooth structure
 - Acidic nature creates microetching of surface, providing micromechanical bonding and ion release
 - Chelating bond: ionic and hydrogen bonds directly to enamel and dentin
- Sometimes uses dentin (cavity) conditioner
 - ~20% polyacrylic acid, removes dentin smear layer
 - Weaker bond than traditional bonding systems, but sufficient




Mechanical properties

- Microhardness increases over time
- Low tensile strength compared to other materials
 - Do not place on high tensile load areas – cusps, incisal edges, etc.



Fluoride release mechanism

- Fluoride within GIC migrates through the water within – can use fluoride from the external environment (saliva)
- Initial release (~4 weeks bolus) is strongest
- Fluoride is concentrated most within the restoration, especially in acidic environments
 - Acid releases more F- because of acid-base reaction in GIC and is concentrated within the restoration
 - Concentration at GIC surface is the least since F- is washed away by saliva



Resin-Modified GICs

- Components
 - Includes water-compatible resin polymers for better handling and mechanical properties
 - Usually hydroxyethyl methacrylate (HEMA) polymer
 - Some safety debates, but no strong evidence against its use
 - Incorporates filler, just like composite, for strength
 - Photoactivator system like composite to assist in handling properties
- Setting reaction
 - Light polymerization occurs first, then the aluminum cross-linking
- Advantages over normal GIC
 - Longer working time with faster set
 - Resin increases strength, resistance to acid



Clinical application of (RM)GICs

- Luting and bonding cements
- Bases or liners in preparations
- Restorations in lower force areas
- Sealants



Part 2 - Cavity lining and base materials

- Liners (<0.5mm thickness) and bases (>0.5mm thickness) are used to protect the pulp from:
 - Thermal stimuli
 - Chemical stimuli
 - Acidity from materials such as zinc-phosphate or zinc-polycarboxylate cements
 - Biological stimuli
 - Microleakage and bacterial entry into crevices
- Material categories
 - CaOH cements
 - Zn-O based cements
 - GI cements
 - RMGICs
 - Visible light-cured resins




CaOH cements (including zinc-oxides)

- Highly basic
- Relatively low compressive strength (~20 MPa), thus mainly used as a liner
- Causes formation of reparative dentin – a calcified layer up to 1.5mm thick
- Consists of 2 pastes:
 - 1: CaOH, zinc oxide, sulfonamide
 - 2: Butylene glycol disilcyate, with other minor chemicals
- Main reaction is a chelating reaction between zinc oxide and butylene glycol disalicylate



Zinc-oxide eugenol cements

- Acid (water) - base (powder) reaction
 - Majority powder, so final product consists of unreacted ZOE held together with salt matrix
- Unmodified ZOE
 - Known for eugenol release - antibacterial properties
 - Poor mechanical strength (compressive strength ~15 MPa) and disintegrates easily – very few uses
- Modified ZOE
 - Additions of hydrogenated rosin to powder and polystyrene / methyl methacrylate to liquid
 - Improves mechanical strength and dissolution resistance
 - Can be used as base or liner (compressive strength ~40 MPa)
- EBA cements
 - ZOE with ethoxybenzoic acid (EBA) - improved mechanical properties to be used as base or liner
- Eugenol reacts with free radicals in resins, thus inhibiting polymerization of resin-based materials – be careful of use



(RM)GIC bases and liners

- Naturally, bonds to tooth structure
 - Can be used under amalgams and composites (composites can bond to GICs as well)
 - To bond with composite, etch the GIC alongside the enamel
 - Etching too long will erode the GIC
 - RMGICs do not need an etched surface due to stronger cohesive strength