



OMFS 711

Principals of Oral Surgery

Dr. Karim El Nokrasky
Assistant Professor
California Northstate University
College of Dental Medicine



Armamentarium for Basic Oral Surgery

Vahid Farahyar, DDS
Assistant Professor
California Northstate University
College of Dental Medicine

4/14/24

California Northstate University College of Dental Medicine



Learning Objectives & Competency Statements

Learning objectives:

- List commonly utilized instruments in basic oral surgical procedures
- Describe different parts of the instruments and instrument handling
- Understand the design of instruments and how they adapt to tooth anatomy



Required Reading

Chapter 7

Instrumentation for Basic Oral Surgery

Hupp, James, R. et al. Contemporary Oral and Maxillofacial Surgery. 7th Edition

Instruments to Incise Tissue

- Scalpel

- Primary instrument for making incisions
- Is composed of
 - A handle
 - A sharp sterile blade
- Single use (plastic handle)
- Reusable handle
 - No. 3 handle
 - most common in oral surgery



Instruments to Incise Tissue

- Scalpel blades


- #15 -most often used blade in intraoral surgeries
 - blade is small and is used to make sulcular incisions and through soft tissue
- #12 is used for
 - Mucogingival procedures
 - Incision on distal of teeth or in maxillary tuberosity area
- #11 has sharp-pointed tip
 - commonly used for small stab incisions (i.e., incising abscess to allow drainage)



Instruments to Incise Tissue

- Inserting/removing scalpel blades

- blade is held along the unsharpened edge, where it is reinforced with a small rib
- male portion of the handle is kept upward
- blade is then slowly slid onto the handle along the grooves in the male portion until it clicks into position



Instruments to Incise Tissue

- *Holding Scalpel*
- Handle is held in pen grasp
- Mobile tissue should be held firmly in place under some tension to maximize cutting efficiency
- For full-thickness mucoperiosteal incision, blade should be pressed down firmly so that the mucosa and periosteum is incised with a single stroke



Instruments to Elevate Mucoperiosteum




- No. 9 Molt periosteal elevator* is the most common instrument
- Pointed end is used to reflect periosteum from teeth, starting at dental papillae
- Broad, rounded end is used to elevate periosteum from bone

Instruments to Retract Soft Tissue


- Retractors

- Provide access and visibility during surgery, and act as shield
- Specifically designed to
 - Retract cheek
 - Retract Tongue
 - Retract mucoperiosteal flap
 - Protect the soft tissue from sharp cutting or rotary instruments



Instruments to Retract Soft Tissue

- *Retractors*




The Henahan (top) and Seldin (bottom) retractors are broader instruments that provide broader retraction and increased visualization

Mission: To advance the science and art of health through education, service, scholarship and social accountability
These teaching materials are the property of California Northstate University and may not be reproduced without permission

Instruments to Retract Soft Tissue


- *Retractors*

- Weider retractor is a large retractor designed to retract the tongue away from the surgical field
- The serrated surface helps engage the tongue so that it can be held securely




Instruments to Grasp Soft Tissue

- *Soft Tissue Forceps*
 - Useful to grab soft tissue during
 - Incision and flap elevation
 - Suturing
 - Locking type
 - Helpful to stop bleeding



Adson tissue forceps **



Stillies pickup



College pliers


Instruments to Grasp Soft Tissue

- Allis Tissue Forceps
 - locking forceps are necessary during removing larger amounts of tissue or doing biopsies
 - locking handle can be held by assistant to apply necessary tension for proper dissection of the tissue



Instruments to Control Hemorrhage

- Hemostats
 - used when pressure alone does not stop the bleeding from a large artery or vein
 - come in a variety of shapes
 - Small or large
 - Straight or curved



Instruments to Remove Bone

- Rongeurs*
 - Bone is cut when sharp blades are squeezed together
 - Two major designs for rongeur forceps
 - Side-cutting forceps
 - Side- and end-cutting forcep



Instruments to Remove Bone

- *Burr and Handpiece*

- Used to remove bone in surgical extraction
- No. 557 or No. 703 fissure burr and No. 8 round burr are commonly used
- For large amount of bone removal, a large-bone burr that resembles an acrylic burr is typically used

557



703



Instruments to Remove Bone

- *Mallet and Chisel*
 - Used for removing lingual tori
 - Sectioning teeth



Instruments to Remove Bone

- **Bone File**
 - Double ended instrument
 - Large end + small end
 - Used for final smoothing of sharp edges or bone spicules before the completion of surgery
 - Teeth of file works well only in **pulling stroke**



Instruments to Remove Soft Tissue From Bony Cavities

- Curette



Angled double ended instrument to remove granulomas, small cysts, and debris from bony defect

Instruments to Suture Soft Tissue

- *Needle Holder*


- has a locking handle and a short, blunt beak
- face of the beak is cross-hatched to ensure a positive grip on the needle



Hemostat

Instruments to Suture Soft Tissue



- *Handling Needle Holder*
 - is held by thumb and ring finger in rings.
 - first finger is placed on top of the instrument and second finger is resting on the side of the instrument to allow better control



Instruments to Suture Soft Tissue


- *Suture Needles (in OS)*

- half-circle or three-eighths–circle
- variety of shapes from very small to very large
- tips of suture needles are either
 - tapered (like sewing needles)
 - triangular tips (cutting needle)
- cutting portion of the needle extends one-third the length of the needle



Instruments to Suture Soft Tissue

- *Suture Needle*
 - is held approximately two-thirds of the distance from the tip of the needle



Instruments to Suture Soft Tissue

• Suture Materials

- Classified
 - By diameter [3-0* / (000)]
 - Resorbability
 - Resorbable
 - Gut/catgut/chromic gut*
 - Polyglycolic/polylactic acid
 - Non-resorbable (silk*, nylon, vinyl, Teflon, etc.)
 - Filament type
 - Monofilament
 - Polyfilament



Instruments to Suture Soft Tissue

- *Suture Scissors*


- Have short cutting tips
- Are held in the same fashion as the needle holder



Spencer Scissors




Dean Scissors



Instruments to Suture Soft Tissue


- *Soft Tissue Scissors*

- Iris
- Metzenbaum



Instruments to Hold the Mouth Open

- *Bite blocks*



Instruments to Hold the Mouth Open

- *Molt/Side Action Mouth Prop*


- May be used
 - In sedated patients
 - In cases of mild trismus




Instruments to Remove Fluids

- *Surgical Suction Tips*

- Remove blood, saliva, irrigation solutions
- Have smaller orifice to remove fluids efficiently



Disposable Surgical Aspirating Tips



Fraser suction tip

Instruments to Hold Towels and Drapes in Position


- Towel Clip
 - ends of the towel clip can be sharp or blunt



Irrigating


- Irrigation Syringes

- Removes debris from surgical site
- Removes debris from bur flutes
- Prevents heat build up and bone damage
- Saline or Sterile Water are used as irrigating solutions




Instruments to Extract Teeth

- *Dental Elevators*



Instruments to Extract Teeth

- *Dental Elevator*
 - Crossbar handle can generate large amounts of force and therefore must be used with great caution



Instruments to Extract Teeth

Types of Elevators

- *Straight Elevators*
 - Most frequently utilized
 - Blade is concave on the working side



Instruments to Extract Teeth

Types of Elevators

- Angled Elevators
 - Most frequently utilized in posterior of the mouth
 - Miller (top image) and Potts (lower images)



Instruments to Extract Teeth

Types of Elevators


- *Triangular Elevators*
(CRYER)
 - Comes in pairs: Right and left or “East and West”



Instruments to Extract Teeth

Types of Elevators


- Pick-Type Elevators
 - Crane pick (top)
 - Root-tip pick or the apex elevator
 - teases small root tips from sockets



Instruments to Extract Teeth

Types of Elevators

- Crane pick



Instruments to Extract Teeth


- Periotomes

- Used for "atraumatic" extraction
- Preserves socket anatomy
- Severs periodontal ligaments



Instruments to Extract Teeth

- Extraction Forceps



Instruments to Extract Teeth


- Extraction Forceps
- *Curved handles*



Instruments to Extract Teeth

Extraction Forceps (Maxillary)

- Upper Universal Forceps (No. 150)
 - Most frequently used to remove single-rooted max. teeth



Instruments to Extract Teeth

Extraction Forceps (Maxillary)

- Upper Primary Teeth Universal Forceps (No. 150S)
 - smaller version of 150 which adapt well to all maxillary primary teeth



Instruments to Extract Teeth

Extraction Forceps (Maxillary)

- *No 1 Forceps*


- Fits roots of maxillary incisors and canines



Instruments to Extract Teeth

Extraction Forceps (Maxillary)

- No. 150A forceps
 - Beak adapt better to maxillary premolar
 - Poor adaptation to anterior teeth



Instruments to Extract Teeth

Extraction Forceps (Maxillary)

- No. 53R and 53L
 - Fit the upper molars
 - One side fits into buccal furcation
 - Concave side fits palatal root



Instruments to Extract Teeth

Extraction Forceps (Maxillary)

- *No. 88R and 88L*
 - Fit the upper molars
 - Accentuated pointed beak
 - Useful when crown of upper molar is severely carious



Instruments to Extract Teeth

Extraction Forceps (Maxillary)

- No. 210S forceps
 - Used for extraction of second molars and erupted third molars that have a single conical root



Instruments to Extract Teeth

Extraction Forceps (Maxillary)

- No. 65 forceps
 - Used to remove broken maxillary molar roots
 - Can be used with ant. Roots too



Instruments to Extract Teeth

Extraction Forceps (Mandibular)

- *No. 151 (Universal Lower Forceps)*
 - Used to remove single-rooted lower teeth
 - Resemble 150, except the beaks point downward



Instruments to Extract Teeth

Extraction Forceps (Mandibular)

- *No. 151S (Universal Lower Forceps for primary teeth)*
 - Resemble 151, scaled down to adapt to lower primary teeth



Instruments to Extract Teeth

Extraction Forceps (Mandibular)

- No. 151A
 - Beaks are parallel for better adaptation to mandibular premolars



Instruments to Extract Teeth

Extraction Forceps (Mandibular)


- *English style of vertical-hinge forceps*
 - Used for single-rooted lower teeth
 - While using, great care must be exercised to prevent root fracture



Instruments to Extract Teeth

Extraction Forceps (Mandibular)


- No. 17 molar forceps
- Beaks enter the bifurcation of lower molars
- Adapts well to right and left molars



Instruments to Extract Teeth

Extraction Forceps (Mandibular)

- No. 87 forceps (cowhorn forceps)
- Heavy beaks enter the bifurcation of lower molars
- Use with caution to prevent root fracture and damage to alveolar bone



Basic Extraction Tray

Weider + Austin + Minnesota retractors

suction tip

No. 9 periosteal elevator

periapical curette

curved hemostat

towel clip

college pliers

straight elevators

suture scissors

Surgical Extraction Tray



Biopsy Tray




Postoperative Tray



Holding Maxillary Forceps



Holding Mandibular Forceps




Holding English Style Forceps



A close-up photograph of a dental procedure. Two hands wearing blue nitrile gloves are using a dental elevator tool to work on a patient's teeth. The background is slightly blurred, focusing attention on the hands and the tool. The text 'Dental Education:' is overlaid in white, serif font, and 'ELEVATOR TECHNIQUE' is overlaid in a larger, bold, white, sans-serif font.

Dental Education:

ELEVATOR TECHNIQUE



English Pattern Mead #1
Upper anterior



English Pattern #7
Upper premolars



English Pattern #17
Upper molars, right, serrated beaks




English Pattern #18
Upper molars, left, serrated beaks



English Pattern #51A
Upper root fragments, serrated beaks



4/14/24





Best Positioning
for Exodontia

Thank You for Your Attention!

He looks confused!
How did he pass his
exams?

